





SABIO140 Automated ESR Analyzer is a bench top system which has complete integration in the hematology laboratories by the direct determination of ESR from blood/EDTA samples with its smart and innovative technology combining "ccd cam" and internationally accepted Westergren method.

SABIO140 is completely cost effective, provides patient comfort and environment friendly since it is using same EDTA tubes coming from CBC tests.

40 sample loading positions, suitablity to continious loading for 23 minutes when you start running the instrument, being maintanence free, being environment-friendly, LCD touch screen and data transmission to LIS are just some of the advanced features of the analyzer which provides laboratory confort and productivity.

- Use of the same CBC tubes. No dedicated citrate tube needed
- Westergren method by CCD camera
- Use of smart cards
- No reagents required, environmently friendly
- 40 sample loading positions, up to 120 tests/hour
- Continuous loading of samples untill the first result
- First result available after 24 minutes and each following result within 5 sec.
- 10 inch LCD touch screen
- User friendly software
- Built-in thermal Printer (optional external printer)
- LIS connection
- Optional Barcode reader



RELIABILITY | PRODUCTIVITY | SATISFACTION



USE OF CBC TUBES

No dedicated ESR citrate tube needed.

No reagents required.



EFFICIENT & USER FRIENDLY
Easy to use. 40 sample positions
and continuous loading



FEEL THE EXCELLENCE...



Technical Specifications

Principles : Westergren method by CCD Camera

Sample Loading Positions : 40, Direct loading of Blood/EDTA Samples

Throughput : Up to 120 samples/hour

Time to First Result : First result in 24 minutes

Time to Result : 5 seconds

Minimum Sample Volume : 2 ml

Sample Test Volume : 2 ml-3 ml

Barcode Reader : External (optional)

Internal Printer : Yes

LIS Communication : Yes

Display : Liquid Touch Screen , 800x600 pixel

Operating Environment : Temperature: 15°C-35°C, Humidity: ≤80%

Power Requirement : A.C. 110-230V ≤300 VA 50/60Hz

Dimension : 690mm(L) x 425mm(W) x 315mm(H)

Weight : 28KG



İzzetpaşa Mah. Vefa Poyraz C. No:5/5, 34387 Şişli İstanbul Phone: +90 212 248 52 64 Fax: +90 212 291 28 65 biota@biota.com.tr

www,biota,com,tr



