

PA54

Specific Protein Analyzer



PA 54

Specific Protein Analyzer

Application	Parameters	Sample Type	Linear Range
Inflammation Monitoring	CRP	Whole blood/serum	3-300mg/L
Cardiovascular Diseases	HS-CRP	Whole blood/serum	0.25-20mg/L
	CCP	Serum	5-100IU/L
Rheumatism Monitoring	RF	Serum/plasma	3-160IU/mL
	ASO	Serum/plasma	20-800IU/mL
Diabetes Monitoring	HbA1c	Whole blood	2-15%
Renal Function Monitoring	CYS-C	Serum/plasma	0-8mg/L
	mALB	Urine	10-200mg/L
Clotting Monitoring	D-Dimer	Plasma	0-20mg/L
	IgG	Serum/plasma	2-40g/L
Immunological Diseases	IgA	Serum/plasma	0.4-8g/L
	IgM	Serum/plasma	0.05-5g/L
	C3	Serum/plasma	0.04-4.5g/L
	C4	Serum/plasma	0.015-2g/L

Tests/Kit: 25T, 2×25T, 50T, 100T

