

DRY BATH INCUBATOR


“ISO 9001:2008 & ISO 13485:2003”
CERTIFIED COMPANY
“CE MARKED PRODUCTS”

Model No : DBI - 11

Mains Input	: 230V/50Hz 110V/60Hz A.C. Supply
W x D x H	: 168 (W) x 232 (D) x 93 (H)
Body	: Metal Body
Weight	: 2.3 Kg. Net
Temp. Controller	: Digital Type
Temp. Range	: 37°C To 70°C
Accuracy	: $\pm 0.5^\circ\text{C}$ @ 37°C, $\pm 3\%$ Except 37°C
Heating Block	: 24 Holes of 12.7mm Diameter or As Selected - any one
Timer	: 1 To 999 (Digital Type)




Block - 1




12.7 MM 24 HOLE

Block - 2




10.5 MM 28 HOLE

Block - 3



18.5 MM - 4 HOLE
12.7 MM - 6 HOLE
15.5 MM - 5 HOLE
10.5 MM - 7 HOLE

Block - 4



12.7 X 12.7 MM
20 SQUAR HOLE



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com